

ACP-700

AUTO CHART PROJECTOR



Brilliantly Conceived & Beautifully Crafted

ACP-700 FEATURES

Leading supplier - The first developed in Korea by the up-to-date technology and experienced production skills all over the world.

HRN
F D A

Clear and Precise Projection Chart

Four different types of chart show clear and bright images for your right vision test.



Comfortable Operation




A remote controller and fast chart change make operator comfortable and stable.

ACP-700 SPECIFICATIONS

Chart	31 Pcs
Chart Selection Speed	Average 0.3sec.
Mask	1 Open 3 Vertical lines 3 Horizontal lines 21 Single letters
Filter	Red / Green
Projection Distance	2.0 ~ 7.0m (3m is standard)
Projection Magnification	30 X (at 5m)
Projection Size	ø 252 mm and 330(W) X 225(H) (at 5m)
Tilt Angle	Ball Joint ($\pm 10^\circ$)
Lamp	6V, 20W Halogen
Auto - Off Funtion	5 ~ 10 Minutes after
Power Consumption	80VA
Power Supply	AC 100 ~ 240V, 50 / 60Hz
CP Dimension / Weight	230(W) X 360(D) X 265(H), 6Kg
RC Dimension / Weight	64(W) X 195(D) X 20(H), 180g
Interface	8-Pin connector for simultaneous control of chart, mask and R/G filter by a motorized refractor's control unit.
Standard Accessories	Each of dust cover Polarized metal screen (400 X 350 mm) Polarized glasser Spare fuses and 2 Batteries Spare lamp
Optional Accessories	Wall bracket, Red / Green glasses

Designs and details can be changed without prior notice for its improvements.

Chart Type (A / B / C / M)



Distributed by