

VOLK SuperField NC®



The Volk Super Field NC® has the same characteristics (magnification) as the popular 90D, but offers a wider field of view. In addition, the working distance relative to comparable fundus loupes is pleasantly large. The glass index as well as the double aspheric glass gives you a particularly good imaging performance.

Features & Benefits:

- Patented double aspheric glass optics provide enhanced imaging
- High index glass provides improved stereopsis and image clarity
- 'Super 90D'. Same magnification with a wider field of view
- Highly versatile lens with multiple contact and non contact options
- Increased working distance compared to competitive fundus lenses

Comes with single lens case

- **Viewing angle: 95 ° / 116 °**
- **Magnification: 0.76x**
- **Laser spot: 1.30x**
- **Working distance: 7 mm**

AVAILABLE IN 7 COLORS (STANDARD IN BLACK):

