

CS@



MIZAR
VIDEOCAMARA DIGITAL
5 MPIXEL

MIZAR
DIGITAL CAMERA
5 MPIXEL



Para Lámpara de hendidura SL9800 | SL9900 | SL99K Elite
for Slit Lamps SL9800 | SL9900 | SL99K Elite

VIDEOCAMARA DIGITAL MIZAR 5 MPIXEL

MIZAR se puede utilizar con las lámparas de hendidura SL9800, SL9900 y SL99K Elite con aumentos de 3X, 5X y Zoom. Esta videocámara destaca por su alto rendimiento y se caracteriza por una excelente reproducción del color. Durante el examen biomicroscópico, el software de aplicación Phoenix permite la visualización de los más mínimos detalles la adquisición de imágenes y vídeos de alta resolución, su gestión y almacenamiento para su posterior revisión.

DIGITAL CAMERA MIZAR 5 MPIXEL

MIZAR can be applied to SL9800, SL9900 and SL99K Elite slit lamps, with 3X, 5X and Zoom microscope. This video camera stands out for its high performance and is characterized by excellent color rendering. While the biomicroscopic examination, the Phoenix application software, allows the display of the smallest details and the acquisition of high-resolution images and videos, their storage and management for further review.



CARACTERISTICAS TECNICAS		TECHNICAL DATA	
Resolución	5Mpx	Resolution	5Mpx
Píxeles de resolución	2448x2048	Pixel resolution	2448x2048
Interfaz	USB 3.0	Camera interface	USB 3.0
Captura de imágenes	35 fps	Frame Rate	35 fps

CPU: I3 o superior (sugerido I5), RAM: 4 GB o superior (sugerido 8 GB), SO OS: Win10x64, Win8x64, Win7x32, Win7x64, GPU: 1 Gbyte no compartido, RESOLUCIÓN: 1280x960 o superior

CPU: I3 or higher (suggested I5), RAM: 4 GB or higher(suggested 8 GB), OS: Win10x64, Win8x64, Win7x32,Win7x64, GPU: 1 Gbyte not shared, RESOLUTION: 1280x960 or higher



CONSTRUZIONE STRUMENTI OPTALMICI
Via degli Stagnacci 12/E
50018 - Scandicci (FI) | Italy
tel: +39 055 722191 | email cs@csitalia.it
WWW.CSOITALIA.IT | WWW.CSOEVENTI.IT

CO104 | Rev. 10/03/2022