

INDIRECT OPHTHALMOSCOPE

Characteristics:

1. The observed images are clear and upward on the fundus.
2. Battery power, convenient for mobile diagnosis.

Specifications:

1. Illumination form: Big spot, Middle spot, Small spot, Red free
2. Diameter of filter: 0.2-0.4MM
3. Illumination brightness: ≥ 80 I X
4. Diopter compensation : -35D~+20D , 24 in total
5. Illumination Source: 2.6V/0.78W, Incandescent Bulb
6. Illumination source : 2# Alkalescent Battery/2# Rechargeable Battery

